

HEALTH REFRIGERATION EQUIPMENT

To generate the best cooling solutions adapted to the specific needs of each client. Technological vanguard, certified quality and the excellence of its products are the guarantees to support its commitment to innovation and that have positioned it as absolute leaders in the health refrigeration equipment

PHARMA

Pharma care

HEALTH

Hospitals
Blood banks

LAB

Laboratories





Quality certificates

- We continue to evolve our commitment to quality, being the first company in the sector to be accredited by a quality management system based on the ISO 9001:2015 standard that encompasses all the machinery of our catalogue.
- The devices described in this catalogue are in accordance with Board Directives DC 2014/30/UE (Electrical Safety) and DC 2014/35/UE (Electromagnetic Compatibility).
- Environmental certificate ISO 14001. "Following its policy of respect and conservation of the environment, Infrico has implemented an environmental management system based on the international standard ISO 14001:2015, through which compliance with all legal requirements is controlled and policies are implemented of reduction of natural resources, recycling and waste management".
- Prevention of risks in the workplace certificate OHSAS: 18001. "Infrico does its work with the maximum guarantees and security; for this purpose, it has implemented a system for the prevention of occupational risks based on the OHSAS 18001:2007 regulation, maintaining an integrated management system: quality-environment-prevention".
- Compliance with the RoHS directive. "We declare that our products comply with the directive ROHS (2011/65/UE), concerning the restriction of the use of certain hazardous substances in electrical and electronic equipment, such as lead, mercury, cadmium, hexavalent chromium, biphenyls and polychromed diphenylethylene in certain concentrations".
- Compliance with regulations on packaging treatment. "We declare that all our packaging is manufactured following international standards for phytosanitary measures, specifically the NIMF 15 standard, which regulates the wooden packaging used in international trade. On all our pallets you will find the seal: Accredits our supplier for the manufacture of packaging treated with thermal shock".
- Compliance with HACCP regulations. "All our equipment can incorporate a HACCP module for compliance with the regulation [CE] no852 / 2004". [optional].
- Compliance with ETL regulations. "On request, Infrico products certified in ETL electrical safety and ETL SANITATION in the United States are manufactured".
- GOST certificate. All our equipment complies with the Russian approval.

PHARMA

Pharma care



According to DIN 58345 all devices must have the following characteristics:

- Operación Temperature Range: $+2^{\circ}\text{C} < T^{\text{a}} < +8^{\circ}\text{C}$
- Enviromental Temperature Range : $+10^{\circ}\text{C} < T^{\text{a}} < +35^{\circ}\text{C}$
- Easy to clean and disinfect.
- Waterproof camera with anti-corrosion protection.
- Mechanical load capacity 100 kg/m^3 .
- Doors equipped with locks.
- Acoustic emission below $<60 \text{ db}$.
- Alarm in case of power failure for a minimum of 12 h .
- Optical - acoustic alarm in case of temperature deviation.
- Potential-free alarm contact.
- Safety element to prevent temperature drop below $+ 2^{\circ}\text{C}$.

HEALTH



According to DIN 58371 all devices must have the following characteristics:

- Operación Temperature Range: $+2^{\circ}\text{C} < T^{\text{a}} < +6^{\circ}\text{C}$
- Enviromental Temperature Range : $+10^{\circ}\text{C} < T^{\text{a}} < +35^{\circ}\text{C}$
- Monitoring of the lower temperature $+ 4^{\circ}\text{C}$.
- Easy to clean and disinfect.
- Waterproof camera with anti-corrosion protection.
- Mechanical load capacity 100 kg/m^3
- Doors equipped with locks.
- Acoustic emission below $<60 \text{ db}$
- Alarm in case of power failure for a minimum of 12 h
- Optical - acoustic alarm in case of temperature deviation $+2^{\circ}\text{C} < T^{\text{a}} < +6^{\circ}\text{C}$
- Potential-free alarm contact
- Safety device to prevent the blood bag from cooling below $+ 2^{\circ}\text{C}$.

LAB



Our range of specific refrigerators and freezers for the different types of laboratories are dimensioned to maintain a homogeneous temperature within the chamber of:

- Refrigeration: $+ 2^{\circ}\text{C} < T^{\text{a}} < + 6^{\circ}\text{C}$

- Freezing: $-15^{\circ}\text{C} < T^{\text{a}} < -86^{\circ}\text{C}$

We take care of:

- Antibodies.
- Patterns.
- Viruses.
- Bacteria.
- Cells.
- Cytokine.
- Enzymes.
- Proteins.
- Biological Products.
- Blood products.
- Reagent test kit
- Cells cultures.
- Clinical samples.
- Reactive.

